

Powder Free Nitrile Examination Gloves BLUE

Item# EM188

Excel's EM188 is a 100% nitrile, latex-free, powder-free, exam glove with a much higher resistance to chemicals than latex. Possesses superior strength and durability with textured fingertips and thickness of 4-5 mil. EM188 is chemo approved.

- 4-5 mil
- Textured Fingertips
- Non-Sterile
- Single-Use Only
- Ambidextrous
- Color Blue
- Barrier Protection
- Glove color contrasts against blood



10 boxes per case * packaged by weight

MADE IN MALAYSIA Tested For Use with Fentanyl Citrate

Item Number	Small	Medium	Large	Extra Large	2XLarge
EM188	EM188-6	EM188-7	EM188-8	EM188-9	EM188-10

Tested for use with Chemotherapy Drugs using ASTM D6978				
Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time			
Carmustine (3.3 mg/ml) Cisplatin (1.0 mg/ml) Cyclophosphamide (Cytoxan)(20.0 mg/ml) Dacarbazine (DTIC) (10.0 mg/ml) Doxorubicin Hydrochloride (2.0 mg/ml) Etoposide (20.0 mg/ml) Fluorouracil (50.0 mg/ml) Methotrexate (25.0 mg/ml) Mitomycin C (0.5 mg/ml) Paclitaxel (Taxol)(6.0 mg/ml) Thiotepa (10.0 mg/ml) Vincristine Sulfate (1.0 mg/ml)	23.3 Minutes, Not Recommended >240 >240 >240 >240 >240 >240 >240 >240			
Fentanyl Citrate and Concentration	Minimum Breakthrough Detection Time in Minutes			
Fentanyl Citrate Injection (100 mcg/2ml)	>240			

Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time			
5-Azacytidine (25.0 mg/ml)	>240			
Carboplatin (10.0 mg/ml)	>240			
Doxetacel (10.0 mg/ml)	>240			
Epirubicin (2.0 mg/ml)	>240			
Gemcitabine (38.0 mg/ml)	>240			
Ifosfamide (50.0 mg/ml)	>240			
Irinotecan (20.0 mg/ml)	>240			
Mitoxantrone (2.0 mg/ml)	>240			
Oncovin (1.0 mg/ml)	>240			
Oxalipatin (5.0 mg/ml)	>240			
Vinorelbine (10 mg/ml)	>240			
Warning: Not recommended for use with Carmustine and Thiotena				

Warning: Not recommended for use with Carmustine and Thiotepa

Gloves used for protection against chemotherapy drug and Fentanyl Citrate exposure should be selected specifically for the type of drug being used. Users are advised to review material safety data sheets to determine the required level of protection. Chemotherapy drug resistant data is available upon request.

Components used in making nitrile glove may cause allergic reactions in some people. If allergic reaction occurs, discontinue use and consult with a physician.

Quality, Consistency, Reliability **Priced Right**

Order Online Today!

Excel Supply Company www.ExcelSupplyCompany.com 2802 Emerald St Milton, WA 98354

Email: Sales@excelsupplycompany.com Ph: 253-896-1195 F: 253-896-1208