Needles and Syringes

High performance products for superior patient care.





Safety and standard pen needles

Our safety and standard pen needles feature thin wall technology for easy flow of medication during injection. The Passive Safety Mechanism on our safety pen needles covers the needle before, during, and after injection and requires no additional steps for activation. They're easy to attach and compatible with most pen injectors available.

Medline Safety Pen Needles

- Safety mechanism covers patient-end needle, reducing chance of needlestick when attached to the pen
- The pen needle is trigger activated and features a red indicator to denote use
- Once activated, the safety mechanism covers the needle on the back end—reducing the chance of needlestick when removed from the pen
- Compatible with most insulin pens including Sanofi Aventis, Lilly, Novo Nordisk, Owen Mumford and Berline Chemie

Ordering information*

Medline's Safety Pen Needle

Item No.	Size	Color	Pkg.
SPN301005	30G x 5 mm	Purple	100/bx, 1200/cs
SPN311005	31G x 5 mm	Purple	100/bx, 1200/cs
SPN311006	31G x 6 mm	Blue	100/bx, 1200/cs
SPN301008	30G x 8 mm	Green	100/bx, 1200/cs

Standard Pen Needles

Item No.	Size	Pkg.
MPHPN2912	29g x 12 mm	100/bx, 10 bx/cs
MPHPN315	31g x 5 mm	100/bx, 10 bx/cs
MPHPN318	31g x 8 mm	100/bx, 10 bx/cs
MPHPN324	31g x 4 mm	100/bx, 10 bx/cs

^{*}For additional sizes and packaging units of measure, contact your Medline Representative.









Standard hypodermic products

Our standard hypodermic needles come in a variety of sizes and styles to meet all your needs.

Ordering information*

Needles

Item No.	Size	Pkg
SYR100167	16g x 1.5"	1000/cs, 10 bx/cs
SYR100185	18g x 1.0"	1000/cs, 10 bx/cs
SYR100187	18g x 1.5"	1000/cs, 10 bx/cs
SYR100195	19g x 1"	1000/cs, 10 bx/cs
SYR100197	19g x 1.5"	1000/cs, 10 bx/cs
SYR100205	20g x 1"	1000/cs, 10 bx/cs
SYR100207	20g x 1.5"	1000/cs, 10 bx/cs
SYR100215	21g x 1"	1000/cs, 10 bx/cs
SYR100225	22g x 1"	1000/cs, 10 bx/cs
SYR100227	22g x 1.5"	1000/cs, 10 bx/cs
SYR100235	23g x 1"	1000/cs, 10 bx/cs
SYR100255	25g x 1"	1000/cs, 10 bx/cs
SYR100257	25g x 1.5"	1000/cs, 10 bx/cs
SYR100253	25g x 5/8"	1000/cs, 10 bx/cs
SYR100302	30g x 0.5"	1000/cs, 10 bx/cs


Blunt Tip Needles

SYR110021	18g x 1"	1000/cs, 100/bx
SYR110022	18g x 1.5"	1000/cs, 100/bx
SYR110021F	Filter 18g x 1"	1000/cs, 100/bx
SYR110022F	Filter 18g x 1.5"	1000/cs, 100/bx

^{*}For additional sizes and packaging units of measure, contact your Medline Representative.

Standard needle color coding

	Color	Size
	Light Yellow	30g
	Orange	25g
	Blue	23g
	Black	22g
	Teal	21g
	Yellow	20g
1	Brown	19g
	Pink	18g
	White	16g



Blunt tip needlesStandard blunt needles
with and without filter

Safety hypodermic products

Nearly one-third of all sharps injuries involve disposable syringes, leading to the transmission of bloodborne viruses that are dangerous to personnel and costly to the healthcare system.¹ The good news is needlesticks are preventable. Medline has variety of easy-to-use safety hypodermic needles to reduce the risk of staff injuries and help improve the safety of your organization.

Needlesticks by the numbers





41%Estimated portion of needlesticks that occur during use¹



1,000Estimated number of sharps injuries per day among hospital personnel¹

Ordering information*

Sure-Snap™ Safety Needles

Item No.	Size	Pkg.
SSN100185	18g x 1"	800/cs, 8 bx/cs
SSN100187	18g x 1.5"	800/cs, 8 bx/cs
SSN100205	20g x 1"	800/cs, 8 bx/cs
SSN100207	20g x 1.5"	800/cs, 8 bx/cs
SSN100215	21g x 1"	800/cs, 8 bx/cs
SSN100217	21g x 1.5"	800/cs, 8 bx/cs
SSN100225	22g x 1"	800/cs, 8 bx/cs
SSN100227	22g x1.5"	800/cs, 8 bx/cs
SSN100235	23g x 1"	800/cs, 8 bx/cs
SSN100237	23g x 1.5"	800/cs, 8 bx/cs
SSN100255	25g x 1"	800/cs, 8 bx/cs
SSN100257	25g x 1.5"	800/cs, 8 bx/cs
SSN100253	25g x 5/8"	800/cs, 8 bx/cs
SSN100273	27g x 5/8"	800/cs, 8 bx/cs
SSN100272	27g x 0.5"	800/cs, 8 bx/cs
SSN100276	27g x 1.25"	800/cs, 8 bx/cs
SSN100302	30g x 0.5"	800/cs, 8 bx/cs

^{*}For additional sizes and packaging units of measure, contact your Medline Representative.






Scan the code to see our Sure-Snap safety needle in action Our Sure-Snap[™]
safety needle and
syringe features flip
technology that engages
a protective shield to
cover the needle after
it's been used.

Safety insulin syringes and needles

Our safety insulin syringe needles can be temporarily covered by pulling the shield forward; after injection, the shield can be permanently locked in place covering the needle. The shield feature can also be used when you're preparing the syringe in one area for safe transportation to the patient in another area.

Ordering information*



For more information or to order, contact your Medline Representative, visit medline.com or call 1-800-MEDLINE.

References: 1. Workbook for Designing, Implementing and Evaluating a Sharps Injury Prevention Program, CDC website. Accessed July 22, 2022. Available at https://www.cdc.gov/sharpssafety/part1TEXTONLY.html#:~:text=Speaker%20Notes%3A%20CDC%20estimates%20that,home%20care%2C%20and%20nursing%20homes 2. Rizk C, Monroe H, Orengo I, Rosen T. Needlestick and Sharps Injuries in Dermatologic Surgery: A Review of Preventative Techniques and Post-exposure Protocols. *J Clin Aesthet Dermatol.* 2016 Oct;9(10):41-49. Epub 2016 Oct 1. PMID: 27847548; PMCID: PMC5104311.

